

## Wheelchair equipment maintenance log

Wheels, brakes, batteries, cushions all need routine checks. *Walk through monthly. Order parts before they fail.*

### MONTHLY CHECK

- Tire pressure / wear
- Brake function
- Battery (if power chair) — full charge cycle test
- Cushion condition — flatness, tears, cleanliness
- Footrests / armrests — secure
- Seatbelt / positioning straps
- Frame for cracks or loose hardware

### SERVICE / REPLACEMENT LOG

DATE	ITEM SERVICED / REPLACED	PROVIDER	COST	NOTES

### VENDOR AND WARRANTY INFO

WHEELCHAIR MANUFACTURER + MODEL

DME PROVIDER NAME + PHONE

WARRANTY INFO / EXPIRATION

*This sheet is for inspiration only. Always consult your medical team for guidance specific to your situation. Inspired Comforts is not a medical provider.*